


### New floor lock system available for Practico operating table

To enhance the stability of the Practico electro-mechanical operating table, a new floor lock option is available. The new floor lock can easily be set by using the brake pedal. With this optional feature, the technical specifications of Practico

changes slightly. The floor lock system can be ordered as a special product, the standard model remains unchanged. Availability from November 2013 onwards.



Enhanced stability for Practico electro-mechanical through the new floor lock system makes the table suitable for even more demanding surgeries.

#### Technical specifications comparison/differences:

	Standard model	Floor lock model
Height/slide	665 - 965 mm	678 - 978 mm
Height/non-slide	595 - 895 mm	608 - 908 mm
Directional castor	Available	Not available
Conductivity	Yes	Yes

The status of the information, specifications and illustrations in this publication is indicated by the date of publication. Whilst the greatest care has been taken over its compilation, Merivaara Corp. is not responsible for any errors or omissions. Merivaara Corp. reserves the right to make changes to the technology, features, specifications and design of the equipment and models without notice. Please check the Merivaara website for the latest updated specifications and other information.

All trademarks are Merivaara property unless otherwise stated. Patents granted or pending.

©Merivaara Corp. 0000EN 112013, printed in markprint