

TRESTON MEDI TREATMENT TROLLEYS

- Easy-clean ABS work surface
- Work surface surrounded by Euro rail for accessories
- Lockable drawers
- Ø 125 mm castors, two of which have brakes. Also available are a 100 mm castor set **LCS100** and 125 mm dual castors **WCS125-meriv**
- The drawers open 80% and can be partitioned with shelf bins and dividers.
- The drawer frame has been powder-coated in light gray.
- The **TS5** label set with 5 pcs 50 x 50 mm labels in five different colors is also available for the drawers:

blue (784 PF Skyblue), yellow (707 PF Primrose Yellow), red (726 PF Medium Red), green (755 PF Light Green) and purple (75 PL Lavender).



1. Treatment trolley NT575D-meriv

- Four lockable drawers
- Interior dimensions of the drawers: 2 pcs 405 x 500 x 110 mm, 2 pcs 405 x 500 x 170 mm
- Exterior dimensions 575 x 640 x 865 mm (W x D x H)
- With threaded rivets
- With threaded rivets and three 120/180/300 mm high drawers **NT575DA-meriv**

2. Treatment trolley NT575

- Three lockable drawers
- Interior dimensions of the drawers: 1 pcs 405 x 500 x 110 mm, 2 pcs 405 x 500 x 170 mm
- Two shelves with pull-out shelf bins
- 6 pcs transparent 6015 shelf bins
- Exterior dimensions 575 x 640 x 995 mm (W x D x H)

3. Treatment trolley NT825D-meriv

- Four lockable drawers
- Interior dimensions of the drawers: 2 pcs 655 x 500 x 110 mm, 2 pcs 655 x 500 x 170 mm
- Exterior dimensions 825 x 640 x 865 mm (W x D x H)
- With threaded rivets
- With threaded rivets and three 120/180/300 mm high drawers **NT825DA-meriv**

4. Treatment trolley NT825

- Three lockable drawers
- Interior dimensions of the drawers: 1 pcs 655 x 500 x 110 mm, 2 pcs 655 x 500 x 170 mm
- Two shelves with pull-out shelf bins
- 10 pcs transparent 6015 shelf bins
- Exterior dimensions 825 x 640 x 995 mm (W x D x H)

TREATMENT TROLLEY ACCESSORIES



1. Upright profile pair NTALX

Code	Dimensions, mm
NTALX	usable height 450

Upright profiles are attached securely to the trolley frame with screws. They provide the base for mounting accessories above the work surface of the trolley.

2. Top shelf NTS

Code	Dimensions, mm
NTS575	425 x 150
NTS825	675 x 150
NTS575-250	425 x 250
NTS825-250	675 x 250

Height-adjustable shelf mounted on upright profiles. The raised edge of the shelf prevents items from dropping off. The corners of the shelf are open to facilitate cleaning.

3. Bin profile NTBP

Code	Dimensions, mm
NTBP575	width 398
NTBP825	width 648

Aluminium rail for mounting stacking and picking bins. Mounted on upright profiles.

4. Euro rail NTE

Code	Dimensions, mm
NTE575	420 x 10 x 25
NTE825	670 x 10 x 25

Aluminium rail for hanging Euro-rail-compatible accessories. Mounted on upright profiles.

5. Euro rail pair for trolley frame NTEF

Code	Dimensions, mm
NTEF	490 x 10 x 25

Additional pair of Euro rails (includes 2 rails) for mounting on the trolley frame.

6. Disinfectant bottle holder DH

Code	Dimensions, mm
DH	

Mounted on a Euro rail. Suitable for Ø 25 mm hand disinfectant bottles. Suitable for 10 x 25 mm rails of any manufacturer.

7. Oxygen bottle holder OCH

Code	Dimensions, mm
UCH	

Available for trolleys with threaded rivets. For 2–3-litre oxygen bottles.

8. Waste bin WB6L and WB10L

Code	Dimensions, mm
WB6L	volume 6 liters
WB10L	volume 10 liters

Attached to the Euro rail from the narrow side of the mount. Includes the mount, bin and lid. Suitable for 10 x 25 mm rails of any manufacturer.

9. Waste bin WB6H and WB10H

Code	Dimensions, mm
WB6H	volume 6 liters
WB10H	volume 10 liters

Attached to the Euro rail from the wider side of the mount. Includes the mount, bin and lid. Suitable for 10 x 25 mm rails of any manufacturer.

10. Bracket for disposable plastic sharps container SCE

Code	Dimensions, mm
SCE23	for 2- and 3-litre bins
SCE47	for 4- and 7-litre bins

Bracket for a sharps container. Mounted on a Euro rail. Fully Sharpsafe® compatible. This product only includes the bracket; disposable sharps containers must be ordered separately. Suitable for 10 x 25 mm rails of any manufacturer.

11. Bracket for examination glove box EGB

Code	Dimensions, mm
EGB	135 x 76 x 221

Bracket for examination glove boxes. Mounted on a Euro rail. Suitable for 10 x 25 mm rails of any manufacturer.

12. Fold-away work surface FDT

Code	Dimensions, mm
FDT	

Available for trolleys with threaded rivets.

13. Drawer divider set

Code	Dimensions, mm
575DIV111	Set for trolley NT575 (3 x 5010-3), gray
825DIV111	Set for trolley NT825 (3 x 5020-3, 1 x 5010-3), gray

Partitioning drawers with shelf bins and dividers keeps items neatly in their places.