

## DORO Adult Horseshoe Head rest

Universal Horseshoe Headrest for adult or pediatric applications in prone or supine positioning. Mounted to a Doro Swivel Adaptor Aluminium or to the Doro Horseshoe Table Adaptor. Gel Pads interchangeable for adult or pediatric application. Permits vertical and lateral adjustments for proper patient positioning.

### Horseshoe system

100 300 900 *Adult headrest with Extension Bar*  
100 300 800 *Pediatric headrest with Extension Bar*



### Pads

100300910 *Adult Right*  
100300911 *Adult Left*  
100300810 *Pediatric Right*  
100300811 *Pediatric Left*

### Cervical Spine Support

100301200

## Contact us

**Merivaara Corp.**  
Tarmontie 2-4  
15860 HOLLOLA  
FINLAND

+358 3 3394 611

[merivaara@merivaara.com](mailto:merivaara@merivaara.com)