

# Orthopaedic extension device 19300

## Especially suitable for treatments when high safe working load is needed

Orthopaedic extension device with high safe working load for Promerix and Practico operating tables for traction of the lower extremities. To be mounted to the foot end of the table.

### Including:

- Extension bars, 2 pieces
- Support feet for extension bars, 2 pieces
- Screw tensioners, 2 pieces
- Traction boots (adults), 2 pieces
- Pelvic support with pad, 1 piece
- Perineal post with pad, 1 piece
- **Includes also leg plates 19348 for supporting legs during the anaesthesia before the final positioning**





## Technical Specifications

|                         |               |
|-------------------------|---------------|
|                         |               |
| Length                  | 820 mm        |
| Height                  | 1857 mm       |
| Height adjustment       | 700 - 1160 mm |
| Safe working load (SWL) | 225 kg        |
| Weight                  | 37.9 kg       |
| Product Code            | 100019300     |

## Contact us

**Merivaara Corp.**  
Tarmontie 2-4  
15860 HOLLOLA  
FINLAND

+358 3 3394 611

[merivaara@merivaara.com](mailto:merivaara@merivaara.com)