


## Special Rotexshoe traction boots

### The Special RotexShoe Traction boots have been adapted to compliment the Merivaara Flexion Extension Device

Stability and flexibility are not always mutually exclusive, but rather properties which are of major importance, particularly in the medical field. They guarantee safe and individual patient treatment, and expand surgical possibilities. The Special RotexShoe Traction boots offer all of that, and easy handling makes your work easier. •

- Increased hygiene standards with the use of disposable padding
- Secure positioning, thanks to a stable, firmly, fixable footrest, padding and belt combination
- Safe grip under heavy loads, tensile strength of up to 80 kg
- No costly additional padding necessary
- Fast positioning time of the extension shoe
- Uniform belt pressure, to reduce shearing forces
- Even pressure distribution on the skin
- Light, highly resistant footrest
- Easy to clean
- Latex free
- Compatible with all standard extension tables
- Easy, intuitive handling



## **Traction boot RotexShoe with triangular interface**

- Order no. RO.0020.2020.MV

## **Disposable inlay incl. mounting brackets**

- 10 pieces
- Order no. RO.0148.2018

## **Traction boot RotexShoe with triangular interface, small**

- Order no. RO.0020.2020.MV.S

## **Disposable inlay incl. mounting brackets, small**

- 12 pieces
- Order no. RO.0148.2018.S

## **Contact us**

**Merivaara Corp.**  
Tarmontie 2-4  
15860 HOLLOLA  
FINLAND

+358 3 3394 611

[merivaara@merivaara.com](mailto:merivaara@merivaara.com)