

DIVIDED LEG SECTION 30248

Gas spring assisted divided leg sections for Promerix and Practico operating tables, with EUR rails.

Section replaces one-piece leg section. Horizontally adjustable between 0–90 degrees, vertically +16 to -90 degrees.



Order no.

100 030 248 Divided leg section 30248

Width	260 mm of one section
Length	732 mm
Length (mattress base)	618 mm
Length (table top)	678 mm
Width (table top)	540 mm
Leg section adjustment	-90 - +16 °
Weight	19.8 kg
Product Code	100060475