

# DIVIDED LEG SECTION 60484

Divided leg section 60484 with solid pelvic section is optimal for Lithotomy cases like laparoscopic procedures. Both foot sections can be adjusted parallel or separately by using table's powered leg joint and gas spring.

Optional leg/arm shield gives more support for legs specially in high lithotomy position.



Order no.

100 060 484 Divided leg section with solid pelvic section, pair

100 060 485 Optional Leg/Arm Shield, single

<b>Width</b>	260 mm of one section
<b>Length</b>	645 mm
<b>Length (table top)</b>	910 mm
<b>Width (table top)</b>	540 mm
<b>Product Code</b>	100060484