


GRAND PROMERIX

Powerful and smooth Grand Promerix is heavy-duty, yet versatile operating table for all kind of surgical procedures. Improved ergonomic design improves the comfort of the whole surgical team. Comfort leads to concentration, confidence and better co-operation. That's efficiency.



C-ARM ACCESS

The renewed slim base gives you perfect C-arm access. Grand Promerix meets even the strictest demands of modern operating theatres easily.



IMPROVED ERGONOMICS

Grand Promerix has lots of improved details for better ergonomics, for example, hidden, smaller castors and inwards cutout for increased foot space. Human-centric design at its best.



INCREASED HYGIENE

Details redesigned with patient safety in mind. Rethought frame make Grand Promerix hygienic and easy to clean.



FLUENT LOCK™

Our Fluent locking system is both user-friendly and reliable. Connect with one click.



HYGIENIC AND EASY TO CLEAN

New bellows and seat section make Grand Promerix hygienic and easy to clean.



NEW DESIGN - IMPROVED DETAILS

Improved details improve hygiene and ergonomics for the operating room staff.

Title	Sliding		Non-Sliding	
Length (table top)	2100 mm		2100 mm	
Width (table top)	540 mm		540 mm	
Width (base)	652 mm		652 mm	
Height adjustment	650 - 1110 mm		610 - 1070 mm	
Max. patient weight	275 kg	600 lb	325 kg	700 lb
Safe working load (SWL)	275 kg	600 lb	325 kg	700 lb
Trendelenburg	-35 °		-35 °	
Reverse Trendelenburg	+35 °		+35 °	
Head rest adjustment	-45 - +45 °		-45 - +45 °	
Back section adjustment	-45 - +80 °		-45 - +80 °	
Leg section adjustment	-105 - +70 °		-95 - +70 °	
Castors size	80		80 mm	
Primary voltage	100 - 240 V		100 - 240 V	
Secondary voltage	24 V		24 V	
Max. power consumption	700 W		700 W	
Frequency	47-63 Hz		47-63 Hz	
Weight	325 kg		300 kg	
Lateral tilt	-25 - +25 °		-20 - +20 °	
Lower frame length	1317 mm		1317 mm	
Max. lifting capacity	490 kg		490 kg	
Product Code	100060000		100060000	