

# MERIMOTE™ WIRELESS REMOTE CONTROL

## Multifunctional, touch-screen remote control for Q-Flow surgical lights and cameras

Merimote™ is a multifunctional remote control for Q-Flow surgical lights and cameras. With its intuitive touch screen Merimote™ guides the user to get full advantage of the working environment.

Operating room staff can control surgical lights and cameras outside the sterile operating area without disturbing the operators.

## Merimote benefits for the user

- Similar easy-to-use user interface as with Q-Flow
- Allows control of Q-Flow and it's camera outside the sterile area
- Wireless controls - no preinstallation required
- Possibility to control several lamp heads in the same operating room in a synchronized way
- Comfortable to use

Available for global deliveries.



Product Code

100060970

---