

# ANAESTHESIA SCREEN III

## Anaesthesia screen III

Adjustable height and angle. Rotates 360 degrees. Made of stainless steel. Universal clamp 10300 is needed for fixing.

## Extension arm for anaesthesia screen III

An extra wing is attached to the anaesthesia screen III, one or two pcs.

Order no.

110111 Anaesthesia screen kit, incl. anaesthesia screen III, 2 x extension arms III and 10300 clamp ( for EUR, US and UK rails)

100 000 III Anaesthesia screen III

100 000 III Extension arm for anaesthesia screen III



Length	590 mm
Width (base)	20 mm
Weight	1 kg
Product Code	110111, 100000111