

ARM AND HAND TABLE 11911

Radiolucent top (350 x 800 mm), light to handle. Swivel-ball leg slides, angles and folds provide unimpeded access. Adjust vertically for tables with X-ray tops. Endrail provided for traction attachment.

Mounts to siderail with two included claps 11043 (EUR/UK), For US rail 2 pcs clamps needs to be ordered separately. Soft pressure management pad improves patient comfort and stability. One support leg, which offers height adjustment from 830 mm to 1320 mm.

Order no.

110911 Arm and hand table 11911 with two clamps 11043
(for EUR, UK or US rails)

100 011 911 Arm and hand table 11911

100 011 043U US-clamp



Width	350 mm
Length	800 mm
Height adjustment	830 - 1320 mm
Max. patient weight	225 kg
Safe working load (SWL)	21 kg
Weight	5.4 kg
Product Code	110911, 100011911