

ARM REST 125

Padded arm board 125 includes a fixing clamp and straps 110.
Adjustable height, gas spring assisted angle adjustment. Rotates
180 degrees. Suitable for EUR-rails and USA-rails.

Optional extra strap:

Straps 110 (pair) for armrests 125, 19180 and 19185.
Weight 0.05 kg, length 800 mm.

Order no.

110125 Arm rest 125 incl. clamp for EUR, US or UK rails
100 000 110 Optional extra strap



| | |
|-------------------------|--------|
| Width | 150 mm |
| Length | 200 mm |
| Height | 640 mm |
| Max. patient weight | 225 kg |
| Safe working load (SWL) | 21 kg |
| Trendelenburg | -36 ° |
| Reverse Trendelenburg | +26 ° |
| Weight | 6.2 kg |
| Product Code | 110125 |