

# ARM REST 19180

Padded arm rest with adjustable height and angle (ball joint). Arm board is easy to operate and adjust by one lever. Rotates 180 degrees. Clamps included.

Optional extra strap:


Straps 110 (pair) for armrests 125, 19180 and 19185.

Suitable for all rails. Weight 0.05 kg, length 800 mm.

## Order no.

119180 Arm rest 19180 with straps 110, incl. clamp (for EUR, US or UK rails)

100 000 110 Optional extra straps (pair)



Armboard positioning

<b>Width</b>	180 mm
<b>Length</b>	600 mm
<b>Max. patient weight</b>	200 kg
<b>Safe working load (SWL)</b>	15 kg
<b>Weight</b>	6.1 kg
<b>Product Code</b>	119180