

CARPA TABLE

Arm and hand table for carpal tunnel and other minor procedures. Unique cantilevered side frames create an extremely rigid surface that is completely radiolucent. Vertically adjustable 50mm. Soft pressure management pad for optimum patient comfort and stability.

Mounts to siderail with two included clamps, 11043 (suitable for EUR and UK rails). For US rail 2 pcs clamps needs to be ordered separately.

Order no:

110912 Carpa table 11912 with two 11043 clamps (for EUR, UK or US rails)

100 011 912 Carpa table

100 011 043U US rail clamp, 1pc



Length	350 mm
Height	650 mm
Max. patient weight	225 kg
Safe working load (SWL)	21 kg
Weight	2.7 kg
Product Code	110912, 100011912