

HAND OPERATING TABLE I 0390

Special shaped radiolucent hand operating table with one support leg. Includes mattress and clamps. Option to use traction device I 0392 and hand fixation device I 0393 (Weinberger-type).

Order no.

I 10390 Hand operating table; Options: EU, US, UK rails

I 00 010 392 Traction device

I 00 010 393 Hand fixation device (Weinberger-type)



Width	850 mm
Height	520 mm
Height adjustment	700 - 1150 mm
Thickness	30 mm
Max. patient weight	200 kg
Safe working load (SWL)	15 kg
Weight	8.8 kg
Product Code	I 10390, I00010390