

HUMERUS TRACTION DEVICE I 10438

I 10438 Humerus Traction Kit, incl. Counter traction post I 1438, Traction unit I 1443, Weinberger hand fixation device I 0393 and two clamps I 0300 for EUR, US and UK rails

I 1438 Counter traction post for humerus, padded (requires a clamp I 0300)

Pad size: length 300 mm, diam. 80 mm.

Total height: 639 mm

Total width: 365 mm

Weight: 2 kg

I 1443 Traction device for intramedullary nailing of the humerus (requires a clamp I 0300).

Total length: 932 mm

Total width: 676 mm. Adjustment range: 0–200 mm

Weight: 3 kg

I 0393 Weinberger hand fixation device, for hand and finger traction. Translucent and sterilizable.

Size: 180 mm x 179 mm

Weight: 0,3 kg



Order no.

I 10438 Humerus traction kit, incl. counter traction post I 1438, traction unit I 1443, Weinberger hand fixation device I 0393, two clamps I 0300 for EUR, US and UK rails

100 011 438 Counteraction post for humerus

100 011 443 Traction device for intramedullary nailing of the humerus

100 010 393 Weinberger hand traction device

100 010 300 Universal clamp

Product Code

I10438, I00011438, I00011443,
I00010393
