

LATERAL ARM BOARD 11802

Arm rest provides arm support for prone, neuro, lateral and chair positions. A wide range of position flexibility, turns 360 degrees. The premium durasoft pad ensures patient comfort. Universal clamp 10300 is needed for fixing.

Order no.

110802 Lateral arm board 11802 with 10300 clamp

100 011 802 Lateral arm board 11802

100 010 300 Universal clamp 10300



Length	150 mm
Height	510 mm
Height adjustment	170 - 480 mm
Max. patient weight	325 kg
Safe working load (SWL)	39 kg
Weight	3.2 kg
Lateral tilt	45 - 45 °
Product Code	110802, 100011802