

ENT-SURGERY HEAD REST 18143

For ENT-surgery. Fits all Merivaara operating tables with "Beach Chair" back rest. For Promerix and Practico the fixing adapter is needed.

Order no.

110143 Head rest 18143 with adapter 60146

100 018 143 Head rest 18143

100 060 146 Adapter for Promerix/Practico



Length	250 mm
Height	110 mm
Trendelenburg	-45 °
Reverse Trendelenburg	+25 °
Weight	3.5 kg
Product Code	110143, 100018143