


ARTHROSCOPIC STRESS POST

A levering fulcrum for improved access to the medical compartment of the knee. Mounts to accessory rail. Universal clamp 10300 is needed for fixing. Pad 210x90 mm. Height adjustment 59–290 mm.

Order no.

I 10530 Arthroscopic stress post incl. clamp 10300 for EUR, US and UK rails

I00 012 530 Arthroscopic stress post



Length	90 mm
Height	470 mm
Height adjustment	59 - 290 mm
Weight	1.5 kg
Product Code	I 10530, I00012530