

# FOOT REST AND TABLE EXTENSION 20114

To prevent the patient from sliding during procedures requiring reverse Trendelenburg. Available with both IS and VEF mattress. Suitable for all Merivaara tables. Universal clamps (2pcs) included.

Mattress size: 276 x 540 x 65 mm.

Order no.

100 020 114 Foot rest/Extension section including mattress

100 000 124 Mattress for the section



Length	276 mm
Length (table top)	270 mm
Height	540 mm
Max. patient weight	225 kg
Safe working load (SWL)	29 kg
Weight	5.1 kg
Product Code	100020114