

# K-TEK™ COMFORT 350 STIRRUPS

The K-Tek Comfort Stirrups for patient positioning for a wide range of gynecologic, urologic and laparoscopic procedures.

The Stirrups increase surgeons' operating room performance by permitting simple adjustment of lithotomy and abduction intra-operatively.

Designed to relieve pressure with soft padding that encapsulates the foot.

K-TEK™ silicone closure straps making cleaning easy.

Flexible finned boots accomodates various leg sizes.

Superior access with lithotomy range of motion from -55° to +85° aids in relieving clinican shoulder/arm strain.

## Order information

I00012806 K-TEK COMFORT 350 STIRRUPS

I00012808 K-TEK™ STIRRUP CART

I00012811 K-TEK™ COMFORT™ 350 STIRRUP BOOT PAD KIT

I00 011 043 RAIL CLAMP, EUR and UK

I00 011 043U RAIL CLAMP, US





<b>Max. patient weight</b>	159 kg	350 lb
<b>Lithotomy Range</b>	+85 - -55 °	
<b>Abduction range</b>	+35 - -9 °	
<b>Product Code</b>	100012806	