

LEG RESTS 20118

Goepel type leg rests are padded and with velcro straps. Wide trajectory thanks to its ball joint. Stainless steel. Universal clamps I0300 (2 pcs) needed for fixing. Fits all Merivaara operating tables.

Order codes:

I10118 Leg rests with straps 20118, pair, incl. 2 x I0300 clamps for EUR, US and UK rails

I00020118 Leg rest 20118 (pair)

I00010300 Universal clamp (1 pc)



Length	320 mm
Height	650 mm
Max. patient weight	225 kg
Safe working load (SWL)	29 kg
Weight	5.1 kg
Product Code	I10118, I00020118