

# STIRRUPS SYSTEM 1258 I

Stirrup system eliminates pressure of the popliteal fossa and perineal nerve. Eliminates uncomfortable rotation, stretching and compression of the knee and hip joints. Rotates 360 degrees. Universal clamp 10300 is needed for fixing.

Order no.

110581 Stirrups with extra soft pads, pair incl. clamps 10300 for EUR, US and UK

100 012 581 Stirrups system, pair

100 010 300 Universal clamp, (1 pc)



Length	265 mm
Height	790 mm
Max. patient weight	225 kg
Safe working load (SWL)	29 kg
Weight	5.4 kg
Product Code	110581, 100012581