

# TOTAL KNEE STABILIZER

Convenient, stable device to support the flexed leg during total knee procedures. Includes pad and clamp I 1043 (EUR/UK). For US rail clamp needs to be ordered separately.

Extra pads I 2586 are available.

Pad I 2586 for total knee stabilizer I 2585. Width 250mm. Weight <0.1 kg.



Order no.

I 10585 Total knee stabilizer I 2585 with I 1043 clamp for EUR, UK or US rails

100 012 585 Total knee stabilizer I 2585

100 011 043U For Us rail clamp

100 012 586 Extra pads

<b>Length</b>	80 mm
<b>Height</b>	250 mm
<b>Height adjustment</b>	55 - 335 mm
<b>Max. patient weight</b>	225 kg
<b>Safe working load (SWL)</b>	53 kg
<b>Weight</b>	1.5 kg
<b>Product Code</b>	I 10585, I00012585