

# CARBON FIBRE BACK SECTION 60502

Merivaara's Carbon fibre sections give surgeons full 360 deg. C-arm imaging during special procedures such as cardiovascular, orthopaedic and trauma operations.

Attached to Promerix operating table's seat section normal or reverse mode.

Includes 80 mm VEF (Visco elastic foam) mattress, which is antistatic. Seams are welded.

## Order information:

I00060502 Carbon fibre back section 60502, incl. VEF mattress, length 1200 mm

I00060503 Mattress for carbon fibre section, 1200 mm, VEF 80 mm



Length	1200 mm
Width (table top)	540 mm
Max. patient weight	275 kg
Product Code	I00060502