

# KNEE SUPPORT 138

Knee support with height adjustable, telescopic support post, radiotranslucent (carbon fibre) and padded, attaches to the adapter. Can be fixed on both sides of the table.

For orthopedic extension device 19300

Weight 2.6kg, adjustable height 105-475mm.



Order no. 100 000 138

<b>Length</b>	500 mm
<b>Height</b>	610 mm
<b>Max. patient weight</b>	225 kg
<b>Safe working load (SWL)</b>	44 kg
<b>Weight</b>	3.6 kg
<b>Product Code</b>	100000138