

ORTHOPEDIC EXTENSION DEVICE I 9300

Suitable for treatment of all fractures in lower extremities.
For Promerix and Practico adapter included.

For Promerix the additional leg section is required. Recommended 40-sections, 60440 Gynecologic or 60240 Standard. Alternatively also 25-sections, 60425 Gynecologic or 60225 Standard can be used.

Adjustments: Extension device: height 700-1160mm, total length 1857mm. Extension bar turns 270 degrees. Pelvic support pad angle: +/-20 degree and +/-40 degree. Perineal post pad rotates 360 degree. Traction device: height 100-350mm, length 165-360mm.



Included in the delivery:

- Table adapter
- Dual joint extension bars (2 pcs)
- Support feet (2 pcs)
- Screw tensioners (2 pcs)
- Traction boots, adult model (a pair)
- Pelvic support with pad
- Perineal post with pad (1 pc)

Order no. 100 019 300

Length	820 mm
Height	1857 mm
Height adjustment	700 - 1160 mm
Safe working load (SWL)	225 kg
Weight	37.9 kg
Product Code	100019300