


# SHORT SUPPORT ARM 19313

For supporting healthy leg with stirrup or leg support.

Accessory holder with clamp for extension bar 19134  
or Accessory rail with clamp for extension bar 19135 needed for  
fixing the additional accessory, e.g. stirrup or leg support.



Order no.  
100 019 313 Short support arm  
100 019 135 Clamp

<b>Length</b>	50 mm
<b>Height</b>	300 mm
<b>Max. patient weight</b>	225 kg
<b>Safe working load (SWL)</b>	58 kg
<b>Weight</b>	1.5 kg
<b>Product Code</b>	100019313