

# LATERAL BRACE I 10652

Lateral brace I 1652 with rectangular pad.

Provide secure and comfortable positioning for hip and other procedures requiring optimum lateral stability. The durasoft pads provide soft cushioning of bony prominences. The swivelling heads conform to anatomical contours, can be securely locked for precise positioning.

Rotates 360 degrees. Pad 220 x 90 mm.

Includes a mounting clamp I 1043. Clamp I 1043 is suitable for EUR and UK rails. If you wish to have a suitable clamp for US rails, please use code E100 011 652 when ordering.



Order no.

100011652 Lateral brace I 1652, incl. 100 011 043 Clamp for EUR and UK rails

E100011652 Lateral brace I 1652, incl. 100 011 043U Clamp for US rail

<b>Length</b>	350 mm
<b>Height</b>	340 mm
<b>Height adjustment</b>	170 - 285 mm
<b>Max. patient weight</b>	181 kg
<b>Weight</b>	1.7 kg
<b>Product Code</b>	100011652