

# SUPPORTING ROLL, HORIZONTAL 10369

Comfortable supporting roll for knee operations.

Total length 454 mm, vertically adjustable 350 mm. Horizontally adjustable 175 mm. Turns horizontally 360 degrees. Pad size: length 300 mm, diam. 133 mm. Requires clamp 10301 for fixing, sold separately.



Order no.

110369 Supporting Roll, horizontal 10369 incl. 10301 clamp

100 010 369 Supporting Roll, horizontal

100 010 301 Clamp 10301

<b>Length</b>	480 mm
<b>Height</b>	440 mm
<b>Weight</b>	3 kg
<b>Product Code</b>	110369, 100010369