

SUPPORTING ROLL, VERTICAL 10368

Total length 662 mm, vertically adjustable 350 mm.
Horizontally turns 360 degrees. Horizontally adjustable 170 mm.
Pad size: length 270 mm, diam. 133 mm.



Order no.

110368 Supporting roll 10368, vertical, incl. 10301 clamp
100 010 368 Supporting Roll, vertical
100 010 301 Clamp 10301

Length	330 mm
Height	660 mm
Safe working load (SWL)	40 kg
Weight	3.2 kg
Product Code	110368, 100010368