

BEDSIDE CABINET ADA

Ada bedside cabinets are modular and durable pieces of patient room furniture.

Hygienic and durable materials and sleek design make the table fit nicely with Merivaara hospital beds. With a wide range of accessories, the Ada cabinet can be modified to meet your needs.

Twin castors make the cabinet easy to move. Castor locks anchor the cabinet firmly in place. Lockable upper drawer helps to keep valuables safe. The optional reading and dinner table is easy to adjust and can be quickly changed from side to side without tools.



- Can be adjusted with a single hand – easy to use

- Stepless height adjustment
- Handedness can be changed easily and quickly without tools
- Individually lockable castors

Title	Bedside cabinet Ada
Depth	600 mm
Height	900 mm
Safe working load (SWL)	50 kg
Castors size	75
Weight	35 kg
Is Lojer Product	I
Product Code	100009160