

TREATMENT TROLLEY NT575, NT575D, NT575D-MERIV AND NT575DA-MERIV

Versatile treatment trolley

Treatment trolleys NT575, NT575D, NT575D-meriv and NT575DA-meriv are particularly suitable for multiple uses.

Trolleys can be used for resuscitation, medicine, sampling and anaesthesia thanks to their wide range of accessories.

Drawers open 80 %. Shelf bins and drawer dividers help to keep things organized. The trolley includes easy-cleanable ABS cover and the work surface is surrounded by Euro rail for accessories.

Depending on the model there is either three or four lockable drawers. Interior dimensions of the drawers are 405 x 500 x 110 mm and 405 x 500 x 170 mm.

Thanks to their modularity and wide range of accessories the treatment trolleys can be used in multiple environments:

- as an anaesthesia trolley in an operating room
- as a nursing and equipment trolley
- as a medicine trolley
- as a sampling trolley
- as a wound care trolley
- as a multipurpose trolley

Multiple accessories can be used to differentiate the trolley according to the needs of the hospital.

Product codes:

NT575: Three lockable drawers and 6 clear shelf bins

NT575D: Four lockable drawers, drawers 2x120 mm and 2x180 mm

NT575D-meriv: Four lockable drawers with threaded rivet, drawers 2x120 mm and 2x180 mm

NT575DA-meriv: Three lockable drawers with threaded rivet, drawers 300 mm, 180 mm and 120 mm



Title	NT575	NT575D-meriv
Width (table top)	575 mm	575 mm
Depth	640 mm	640 mm
Height	995 mm	865 mm
Product Code	NT575, NT575D, NT575D-meriv, NT575DA-meriv	NT575, NT575D, NT575D-meriv, NT575DA-meriv