

HOSPITAL BED CARENA

The ergonomic choice for both patients and nurses

Carena is a revolutionary hospital bed with unparalleled ergonomics for the patient and for the nurse. The needs of nursing staff as well as an increase in the average weight of the population and the relative number of elderly people have been taken into consideration in its design.

Smaller patients will love it – the minimum height makes it easy for them to get in and out of the bed on their own.

Additionally, the hospital bed accommodates also heavier patients the maximum safe working load (SWL) being 330 kg.

This next-generation hospital bed is equally suited for use in day surgery, emergency room, long-term or intensive care or nursing homes.





- High safe working load 330 kg
- Unique lower and upper height 310–890 mm (with 125 mm castors)
- Comfortable and ergonomic for both patient and nurse
- Light weight and easy transportation
- Electric and manual Trendelenburg
- Quick release (CPR) for back section (optional)
- Innovative structure and design, patented technology
- Long service life-cycle and convenient maintenance
- Comply with EN 60601-2-52

Title	Hospital bed Carena
Width	800 / 850 / 900 mm
Length	2154 mm
Height adjustment	315 - 815 mm
Max. patient weight	300 kg
Safe working load (SWL)	330 kg
Trendelenburg	12 °
Reverse Trendelenburg	6 °
Back section adjustment	0 - 70 °
Leg section adjustment	-24 - +29 °
Protection class (IP)	66
Is Lojer Product	I
Product Code	100000360

