

PATIENT TROLLEY EMERGO

Comfort and reliability in acute care and recovery

Merivaara's multipurpose Emergo patient trolleys / stretchers are used for patient transportation, emergency, recovery, ICU and minor procedures. Their safe and sturdy construction makes Merivaara Emergo patient trolleys ideal for heavy use. The versatile range of options and adjustments ensure efficient patient transfer, patient comfort and ergonomic working positions for health-care staff.

The adjustable mattress base provides comfortable positions for patients to recover eliminating unnecessary patient transfer from the trolley (e.g. in sitting position). Large, easy rolling castors provide excellent mobility and the optional fifth wheel improves overall mobility reducing the strain for the staff.

Emergo multipurpose patient trolley has many attachment places for the IV-rod and other accessories. The wide range of accessories available will serve the demanding needs of nursing staff for many years to come.

Designed and manufactured in Finland, Emergo is a long lasting, good-value patient trolley for any clinic or hospital.





- Modular structure, can easily be equipped according to your needs
- Great to manoeuvre even in tight spaces, fifth wheel improves mobility even further
- Light, gas spring-assisted adjustments
- A wide range of accessories available, including ABS storage boxes
- X-ray platforms available for 2 and 3 sectional mattress bases.
- 6370: X-ray platform opening sizes are: side width 49 cm (leg section), 35 cm (middle section), 60 cm (back section) and depth 47 cm (all sections)
- 6270 and 6280: X-ray platform opening sizes are: side width 47,3 cm (legs), 51 cm (middle), 74,7 cm (back section) and depth 43,9 cm (legs), 47,9 cm (middle), 43,5 cm (back section)

Title	6250	6270	6280	6350	6370
Width	750 mm	750 mm	750 mm	750 mm	750 mm
Length	2075 mm	2075 mm	2150 mm	2075 mm	2075 mm
Height adjustment	560 - 930 mm	625 - 995 mm	625 - 995 mm	560 - 930 mm	625 - 995 mm
Safe working load (SWL)	225 kg	225 kg	225 kg	225 kg	225 kg
Trendelenburg	-14 °	-14 °	-14 °	-14 °	-14 °
Reverse Trendelenburg	+7 °	+7 °	+7 °	+7 °	+7 °
Back section adjustment	0 - 70 °	0 - 70 °	0 - 70 °	0 - 70 °	0 - 70 °
Leg section adjustment	0 °	0 °	0 °	0 - 45 °	0 - 45 °
Castors size	150 / 200	150 / 200	150 / 200	150 / 200	150 / 200
Castor type (antistatic, locking)	antistatic	antistatic	antistatic	antistatic	antistatic
Mattress base (sections, pcs)	2	2	2	3	3
Mattress base (width)	705 mm	705 mm	705 mm	705 mm	705 mm
Is Lojer Product	I	I	I	I	I
Product Code	100006250, 100006270, 100006280, 100006350, 100006370	100006250, 100006270, 100006280, 100006350, 100006370	100006250, 100006270, 100006280, 100006350, 100006370	100006250, 100006270, 100006280, 100006350, 100006370	100006250, 100006270, 100006280, 100006350, 100006370