

WASHABLE HOSPITAL BED CARENA

Purely Finnish design

Merivaara's popular Carena bed has all the features required by healthcare professionals, both from the point of view of ergonomics and patient safety. Washable Carena complies with the strictest standards and hygienic criteria.

The new special Carena model is suitable for the demanding conditions imposed by machine washing with water. Additionally, its antibacterial coating reduces bacterial growth, also on the frame. Therefore, Carena is an ideal choice, both in primary care and ER and in long-term care.

Did you know...?

A bed in the lake?



Special features of the machine washable bed

- Corrosion protection of the bed frame prevents rusting
- Ingress protection class IP66W increases the service life of electrical components
- Antibacterial coating restrains bacterial growth
- Furthermore, dual castors improve the bed manoeuvrability, making transportation easier

General features of Carena

- Four-sectional mattress base, electrical adjustments
- Safe working load (SWL) 330 kg – suitable also for large patients
- Height adjustment 340–855 mm (with 150 mm castors)
- Wide variety of positions that support treatment and provide unparalleled patient and nurse ergonomics
- Exceptional manoeuvrability and light brake system
- Wide array of accessories
- Patented features

Title	Washable Carena
Width (base)	850 mm
Height adjustment	340 - 835 mm
Max. patient weight	310 kg
Safe working load (SWL)	330 kg
Trendelenburg	12 °
Reverse Trendelenburg	6 °
Back section adjustment	0 - 70 °
Leg section adjustment	24 - 29 °
Mattress base (sections, pcs)	4
Product Code	100000360