

LEG RESTS, POWERED

Powered leg rests are suitable with Promerix and new, smarter Practico operating tables, replacing a leg section.

Leg rests provides excellent patient comfort and precise surgical positioning. The VEF-pad reduces popliteal pressure. Velcro straps for security. Integrated clamps included. Adjustment ranges: Support arm (vertical) -95 - +70 degrees, supports (horizontal and vertical) 360 degrees.



Order no. 100 060 482

Length	600 mm
Height	345 mm
Max. patient weight	275 kg
Weight	3.0 kg
Is Lojer Product	0
Product Code	100060482

