

ARM REST 151 FOR OPHTHALMOLOGY

Arm rest pair 151 for ophthalmology is designed to support surgeon's arms during operation. Arm board turns 360 degrees. Height adjustment with joint screw on the side of the aluminium clamp. Fits all Merivaara operating tables.

Order no.

100 000 151 Arm rest 151, pair, incl. clamps for fixing



Length	85 mm
Height	300 mm
Height adjustment	120 - 320 mm
Safe working load (SWL)	29 kg
Weight	2.4 kg
Is Lojer Product	0
Product Code	100000151

