

BOWL 6L, FOR PROMERIX

Meant to be used in Promerix table attached under the 40-gynaegologic section for Uro and Gyn procedures.

61581 without drainage, 61582 with drainage. Stainless steel. Holder for 6l bowl 61580. Attached to gyn section 60440 or 60425.

Bowl size: 396 x 265 x 127 mm

Weight: 2–2.5 kg

Drainage: silicon tube diameter 19 mm, length 1000 mm

Order no.

100 061 580 bowl holder

100 061 581 holder + bowl without drainage

100 061 582 holder + bowl with drainage



Is Lojer Product

0

Product Code

100061581, 100061582

