

HEAD REST FOR NEUROSURGERY I 44

Neurosurgical head rest I 44 is padded and fits all Merivaara operating tables. For fixing Promerix and Practico operating tables, adapter 60146 is required.

Order no.

110144 Head rest for neurosurgery I44 with 60146 adapter

100 000 I44 Neurosurgical headrest I44

100 060 I46 Adapter for Promerix/Practico



Height adjustment	0 - 290 mm
Weight	3.4 kg
Is Lojer Product	0
Product Code	110144, 100000144

