

HEAD REST HELMET TYPE

Helmet type head rest 30128 for shoulder and upper arm procedures with "Beach Chair" back rest.

Unique sliding function with double joint mechanism, very wide adjustments in all directions. Fits Promerix 60525 and Practico 30125 back rests.

Rotates 360 degrees, adjustments 180 degrees.
Vertical adjustment: 65 mm.

Order no.

145128 special head rest for shoulder surgery - helmet model
30128 with adapter 60146

100 030 128 special head rest for shoulder surgery - helmet model



Length	450 mm
Height	605 mm
Max. patient weight	225 kg
Safe working load (SWL)	12 kg
Weight	4.6 kg
Is Lojer Product	0
Product Code	145128, 100030128

