

KNEE CRUTCHES 11181

Provides excellent patient comfort and precise surgical positioning. The "Doughty soft" pad reduces popliteal pressure. Velcro straps for security. Universal clamp 10300 is needed for fixing. Knee crutches turn 45 degrees and rotate 360 degrees.



Order no.

110181 Knee crutches, pair, with extra soft pad and straps 11181 with 2 x 10300 clamps for all rails

100 011 181 Knee crutches 11181 (pair)

100 010 300 Universal clamp

Length	320 mm
Height	650 mm
Height adjustment	450 - 520 mm
Max. patient weight	225 kg
Safe working load (SWL)	29 kg
Weight	2.8 kg
Is Lojer Product	0
Product Code	110181 , 10001181

