

# LIGHT LEG SECTION FOR MOTORIZED LEGS

Light leg section improves work ergonomics of the OR staff. It is very light and easy to handle. Light weight leg section attaches to the operating table with easy locking mechanism.

Light Leg Section for motorized legs, 100060204, is suitable to be used with Promerix operating table and New, smarter Practico with motorized leg joints. Includes 80 mm thick, antistatic VEF (viscoelastic foam) mattress.

Order code:

100060204 Light Leg Section for motorized legs, VEF mattress

For non-motorized legs joints, see Light Leg Section for non-motorized legs >>



Length	630 mm
Width (table top)	540 mm
Height	76 mm
Safe working load (SWL)	40 kg
Weight	4.9 kg
Is Lojer Product	0
Product Code	100060204

