

OPERATING TABLE PRACTICO MANUAL

For a wide variety of surgical procedures

The new hydraulic Practico is a versatile manually operated table model that enables a wide range of configurations for general and special surgery.

The table accommodates also larger patients, and offers excellent access to the operated site thanks to its design and adjustment options.

Additionally, the rust proof table top design makes the table suitable also for demanding climate conditions.



Practico manual features

- Excellent access for the surgical team to the patient
- Versatile – use the same table for a wide variety of surgeries
- Seat section is equipped with GYN cutout
- Faster turnaround times thanks to easy set-up
- Compatible with many of the same accessories and table top sections as other Merivaara tables (Practico electro-mechanical and Promerix)
- C-arm and X-ray compatible
- Suitable for humid weather conditions thanks to the corrosion-free material

Width	620 mm
Length	1950 mm
Width (table top)	594 mm
Height adjustment	710 - 1020 mm
Safe working load (SWL)	180 kg
Trendelenburg	25 °
Reverse Trendelenburg	18 °
Head rest adjustment	-45 - +45 °
Back section adjustment	-30 - +70 °
Leg section adjustment	-90 - +15 °
Mattress base (width)	540 mm
Weight	135 kg
Lateral tilt	-20 - +20 °
Max. lifting capacity	400 kg
Product Code	100030010