

PEDIATRIC SECTION 60506

Pediatric section is designed for children's surgery. The table top is narrower than other sections enabling surgeons and health care staff to get closer to the patient.

Max patient weight is 40 kg. Includes 80 mm thick antistatic and viscoelastic mattress.

The section is meant only for Promerix operating table's seat section's joints on back or leg side.

Order nr. I00060506



Title	Pediatric section 60506
Length (table top)	1340 mm
Safe working load (SWL)	40 kg
Mattress base (width)	350 - 540 mm
Weight	18 kg
Product Code	I00060506

