

TOTAL KNEE STABILIZER

Convenient, stable device to support the flexed leg during total knee procedures. Includes pad and clamp I 1043 (EUR/UK). For US rail clamp needs to be ordered separately.

Extra pads I2586 are available.

Pad I2586 for total knee stabilizer I2585. Width 250mm. Weight <0.1 kg.



Order no.

I 10585 Total knee stabilizer I2585 with I 1043 clamp for EUR, UK or US rails

100 012 585 Total knee stabilizer I2585

100 011 043U For Us rail clamp

100 012 586 Extra pads

Length	80 mm
Height	250 mm
Height adjustment	55 - 335 mm
Max. patient weight	225 kg
Safe working load (SWL)	53 kg
Weight	1.5 kg
Is Lojer Product	0
Product Code	I 10585, I00012585

