

MONITOR ARMS

Monitor arms for diverse needs

With Merivaara monitor arms, equipment is easily adjustable to correct viewing and operating position and can be moved away if necessary. Soft adjustments and stable operation ensure that monitor arms remain in the desired position.

One or more displays can be installed on Merivaara monitor arms. All monitor arms are available as wall or ceiling mounted versions.





Merivaara monitor carrier system benefits

- Soft adjustments and stable operations
- Flexible tilt angles provide optimal monitor positioning of adjacent screens
- Wall- or ceiling-mounted spring arms available
- With power and data cables fully enclosed for superior hygiene
- Multiple configurations available
- One or more displays can be installed
- Tailored solutions available with support from Merivaara professionals

Title	Acrobat 2000	Acrobat 3000	OndaSpace
Maximum load	21 / 2x18 kg duo	30 kg	40 kg
Weight ranges	1–3,5 / 3,5–7 / 3–6 / 7–12 / 8–13 / 9–15 / 12–18 / 13,5–21 kg	7–11 / 8–12 / 12,5–18 / 17–26 / 23–30 kg	8–16 / 15–32 / 20–40 kg
Lead-in cables	Esim. 2 kpl max 9 mm	Esim. 1 kpl 9 mm + 2 kpl 3,5 mm + 1 kpl 6 mm ja 2 kpl 5 mm	Esim. 3 kpl max 12 mm + 2 kpl 6,5 mm
Spring arm length	910 mm	930 mm	988 mm
Optional fixed arm length	600 / 750 / 950 mm	700 / 850 / 1000 / 1150 mm	700 / 850 / 1000 / 1150 mm
Additional information	Yhdelle monitorille, saatavissa myös kahdella monitorivarrella	Yhdelle/kahdelle monitorille	Yhdelle/kahdelle monitorille
Product Code	Monitor carrier systems	Monitor carrier systems	Monitor carrier systems

