

Anaesthetic trolley 650

Anaesthetic trolley 650 is made of stainless steel.

- 6 Drawers, with easy sliding motion, one with lock
- Wire basket for catheters and hoses at the back of the trolley
- Side rail for fixing suction bottles, scopes, blood pressure meter, catheters etc.
- Castors \varnothing 125 mm, antistatic, 2 castors with brakes



Technical Specifications

Width	540 mm
Length	785 mm

Height	850 mm
Safe working load (SWL)	100 kg
Product Code	100000650

Contact us

Merivaara Corp.
Tarmontie 2-4
15860 HOLLOLA
FINLAND

+358 3 3394 611

merivaara@merivaara.com