


## 4K/UHD for wall mounting

Merivaara offers a wide range of medically certified monitors suitable for visualizing the content during medical procedures. We offer also a range of compatible accessories and even customized mounting solutions to meet your demands.



### Available models

**SUPPLIER CODE**

[FM-C5501DV](#)

[FM-D5801DV](#)

[MBT-1000.55.2](#)

[MDAC-8355 SPIB](#)

[MDAC-8355 SPIW](#)

[MDSC-8255 LED](#)

[MDSC-8358 RL](#)

[MDSC-8358 RLG](#)

**PRODUCT CODE**

A43818900

A43819100

A43861700

A44043900

A44044000

A43854600

A43854800

A43854900

## Contact us

**Merivaara Corp.**  
Tarmontie 2-4  
15860 HOLLOLA  
FINLAND

+358 3 3394 611

[merivaara@merivaara.com](mailto:merivaara@merivaara.com)