


Q-Flow compatible Full HD monitors

Merivaara offers a wide range of medically certified monitors suitable for visualizing the content during medical procedures. We offer also a range of compatible accessories and even customized mounting solutions to meet your demands.



Available models

SUPPLIER CODE

FS-L2403D
FS-L2702D

PRODUCT CODE

A45000185
A45000472

Contact us

Merivaara Corp.
Tarmontie 2-4
15860 HOLLOLA
FINLAND

+358 3 3394 611

merivaara@merivaara.com